



User Manual

SwissSTIM Physio V 7.0

2005

I. Introduction

The **Swiss_{STIM}™ Physio** is manufactured/distributed by VALMED SA, Sion, Switzerland.

The **Swiss_{STIM} Physio** is manufactured in accordance with the requirements of European Safety Standards **EN 60601-1**, **EN 60601-2-10** and meets all requirements of the American Safety Standards for Transcutaneous Stimulators as set forth in **ANSI/AAMI NS4 – 1985**. It is approved by the United States Food and Drug Administration (FDA) (K022175).

The **Swiss_{STIM} Physio** is a Class II Medical Device and conforms to the requirements of European Directive CEE 93/42 and holds certificate number **CE 0476**.



BF device IAW IEC 601-1 "External Placement on Body".



Physio . Read this *User Manual*, especially SECTIONS 3 and 4 **before** using the **Swiss_{STIM}**

Manufacturer is not liable for any damage caused by the improper use of the **Swiss_{STIM} Physio** neuromuscular stimulator.

II. Table of Contents

I.	Introduction	2
II.	Table of Contents	3
III.	Indications, Contraindications, Warnings, Precautions and Guidelines	4
IV.	Terminology, Safety, Treatment and Electrode Data	6
V.	Swiss_{STIM} Physio Unit Components	9
VI.	Operating Instructions	10
VII.	Program Descriptions	15
VIII.	Useful Information	20
IX.	Swiss_{STIM} Physio Technical and Safety Data	21
X.	Electrode Placements and Program Guide	23
XI.	Program Reference Guide	26

III. Indications, Contra-Indications, Warnings, Precautions and Guidelines

A. Indications

The **Swiss_{STIM} Physio** is intended for use by licensed health practitioners, patients under the care of licensed health professionals and is limited to prescription sales. The **Swiss_{STIM} Physio** treatment programs are designed and intended for stimulation of all parts of the body EXCEPT transthoracically, the head and the front part of the neck.

The specific indications for use include:

- Prevention or retardation of muscle disuse atrophy
- Relaxation of muscle spasms
- Increasing local blood circulation
- Muscle re-education and strengthening
- Immediate post-surgical stimulation of calf muscles to prevent venous thrombosis
- Maintaining or improving range of motion of extremities

B. *Medical Contraindications*

Do **NOT** use the **Swiss_{STIM} Physio** on patients who have:

- An implanted cardiac pacemaker, defibrillator or other implanted electronic or metallic devices
- Any type of cancer in active metastatic phase.

C. *Warnings for Medical Professionals*

- Do not use in the presence of functioning high frequency electrosurgery devices
- Do not use in the immediate vicinity (< 0.5 meter/2 feet) of active microwave ovens
- Long-term effects of chronic electrical stimulation are not known.
- Keep out of the reach of children
- Apply electrodes ONLY to normal, intact, clean skin. Do not apply electrodes over open wounds or over swollen, infected or inflamed areas or skin eruptions, e.g., phlebitis, thrombophlebitis, varicose veins, etc.
- **DO NOT APPLY STIMULATION:**
 - Over frontal area of the neck (near carotid sinus nerves).
 - Over the neck or mouth. Severe spasms of the laryngeal and pharyngeal muscles may occur with contractions strong enough to close the airway and/or cause difficulty in breathing. Stimulation over the neck could also have a adverse effects on the heart rhythm or blood pressure.
 - Transcerebrally or transthoracically
 - Over, or in proximity to, cancerous lesions.
 - On patients who have a hernia (abdominal or inguinal)
 - To the frontal, laryngeal and temporal regions of the neck.
- Advise patients to never use the **Swiss_{STIM} Physio** while driving, operating machinery or during activities in which involuntary muscle contractions may endanger the user or others.
- The effects of stimulation of the brain are unknown. Therefore, do not apply stimulation across the head and do not place electrodes on opposite sides of the head.
- Do not allow use of the **Swiss_{STIM} Physio** in the bath or shower.
- Patients with suspected heart problems or epilepsy must obtain appropriate medical advice.
- Never allow use of the **Swiss_{STIM} Physio** while sleeping.
- NEVER immerse the **Swiss_{STIM} Physio** unit in any liquid.

D. Health Precautions

- The safety of neuromuscular stimulation during pregnancy has not been established.
- Use caution when/if:
 - Patient has skin areas that lack normal sensation.
 - Following surgical procedures if muscle contractions might impede the healing process.
 - Over a menstruating or pregnant uterus.
 - There is a tendency to bleed internally following acute trauma or fracture
- Place electrodes in accordance with illustrations in this User Manual.
- Ensure that patient extremities are isometrically fixed (braced) during treatment sessions to prevent movement that results from stimulation.
- This unit should not be used while driving, operating machinery or during any activity in which involuntary muscle contractions may place the user at undue risk of injury.
- Some users may experience skin irritation or hypersensitivity due to the electrical stimulation or the conductive medium.
- Application of moderate heat (thermal wrap) to muscles as well as moistening skin prior to treatment improves treatment efficacy; use of cold packs on treated muscles after treatment is likewise recommended.
- This unit should only be used with the leads, electrodes and accessories provided by the manufacturer.

E. DO NOT use the Swiss_{STIM} Physio on patients who

- Have an implanted cardiac pacemaker, defibrillator or other implanted electronic or metallic devices
- Are pregnant
- Have cardiac problems or cardiac disease
- Have epilepsy
- Have abdominal or inguinal hernia

F. Electrode Guidelines

The PalsFlex electrodes that are supplied as a standard accessory with the **Swiss_{STIM} Physio** mold easily to body contours and are reusable. The specific instructions for electrode use are indicated on the factory sealed storage pouch. The recommended sizes are oval 3" by 5" electrodes for large areas (e.g., leg muscles) and round 3" electrodes for smaller areas such as forearm muscles.

- Use only skin pads (electrodes) supplied by manufacturer. Other electrodes may not stimulate correctly and may render the treatment ineffective.
- Apply electrodes only to clean, intact, normal skin.
- Do not apply electrodes over open wounds, inflamed, swollen or infected skin area or over any skin eruptions such as varicose veins, phlebitis, etc.
- Do not share electrodes with other patients. Multiple patient use may result in adverse skin reactions.
- Replace self-adhesive electrodes when they do not adhere (stick) firmly to the skin

G. Potential Adverse Reactions

Skin irritation and burns beneath the electrodes has been reported from use of some neuromuscular stimulators.

IV. Terminology, Safety, Treatment and Electrode Data

A. Terminology

Human physical movements are determined by muscle actions on the skeletal system. A muscle is shortened (contracts) when it is “voluntarily” activated by the brain. This shortened muscle exerts force on the attached bone(s), causing movement with skeletal joints acting as pivots. All muscle contractions are controlled by the brain through very efficient electrical stimulation signals that travel through motor nerves. A signal will activate a discrete number of muscle cells known as a “**motor unit**”; a motor unit is a single motor nerve and the muscle cells connected to it. Groups of motor nerves and muscle cells are known as “**motor points**”.

The full contraction of a muscle typically involves multiple motor units acting simultaneously; the force of contraction is proportional to the number of motor units activated. The gradual activation of motor units enables smooth and controlled development of force; this process is known as “**spatial summation**”.

Neuromuscular stimulation using the **Swiss_{STIM} Physio** achieves similar muscle movement without use of brain signals. These movements are termed “involuntary”.

When a single electrical impulse of adequate intensity is applied **externally (involuntarily)** to a group of neuromuscular junctions (**motor points**) in a muscle, the result is a single short contraction (twitch) in that muscle. When these single twitches are repetitive and the repetition rate exceeds 10 twitches per second, the contractile force of each succeeding twitch adds to the preceding twitch an additional degree of contraction, resulting in a higher overall force and, hence, muscle contraction.

This “effect” is termed “**temporal summation**”.

The minimal repetition rate at which twitch contractions “fuse” together is termed the “**tetanzation frequency**”, typically in the range of 25 to 50 impulses per second, depending upon the specific muscle. Tetanzation frequencies are used to create electrically stimulated contractions known as tetanic muscle contractions.

During natural (voluntary) muscle contractions, the force is partly the result of summation of single repetitive twitches (temporal summation) but also a function of the total number of motor units that have been activated by brain signals. Therefore, a voluntary **muscle contraction force** is a result of both temporal and spatial summations.

During neuromuscular electrical stimulation, the muscle contraction force is highest at the tetanic stimulation frequency due to temporal summation. Contraction force also depends upon the total number of motor units activated, which, in turn, are dependent on placement of electrodes and their distance to the motor points on a muscle.

B. Safety

The **Swiss_{STIM} Physio** is designed to provide a totally safe treatment without sacrificing effectiveness. In this respect, the **Swiss_{STIM} Physio** is non-pareil. Further, to ensure safety, durability and efficiency, only premium electrical components are used.

All stimulation impulses generated during **Swiss_{STIM} Physio** treatment sessions carry such minimal amounts of electrical energy that they are unlikely to produce any adverse effects when

the stimulator is used in accordance with this manual. It is important, however, to emphasize that no neuromuscular stimulator, including the **Swiss_{STIM} Physio** should be used by a patients with an implanted cardiac pacemaker and that the safety standards for pregnant women has not been established.

The **Swiss_{STIM} Physio** is designed so that even improper or accidental application of the stimulator will not produce cardiac rhythm disturbances (this does not apply to persons who have implanted cardiac pacemakers). This safety factor is due to the minimal electrical charge of the stimulating impulses, which, under all conditions, do not exceed 24 microcoulombs, zero net current. This charge is below the standard of 25 microcoulombs established in the Association for Advancement of Medical Instrumentation (AAMI) for cardiac rhythm disturbance safety (AAMI/ANSI NS-4-1985).

Many, if not most, individuals harbor apprehension, anxiety or even fear regarding electricity. It is therefore very important that patients understand how completely safe treatments with the **Swiss_{STIM} Physio** are.

Educating the patient before treatment is therefore a recommended first step!

The effective value of the stimulation voltage a patient will experience during treatment with the **Swiss_{STIM} Physio** is very low, below 5 volts (root mean square) at the maximum setting of the intensity knobs.

C. Practical Treatment Suggestions

During neuromuscular stimulation, the number of motor units in a muscle that are activated depends upon the stimulus impulse energy. This imparted impulse energy, in turn, is a function of:

- The current intensity selected by the practitioner for the patient;
- Skin and electrode electrical resistance; less resistance results in more delivered energy, and;
- Electrode placement; the closer electrodes are to motor points, the higher the energy.

The intensity of the stimulation should be at a level where the involuntary muscle contractions are visible and felt by the person to be treated. The controls on the **Swiss_{STIM} Physio** are easily adjusted to achieve this. It is likewise true that the stronger the muscle contractions are, the higher the effectiveness of the therapy. This has to be tempered, however, with a person's comfort level. For comfortable and effective stimulation, the delivered energy should be optimized. This can be achieved by a combination of the following steps:

- Ensure that the person's skin is clean;
- Apply heat (for example, a thermal wrap) to the skin to increase local blood circulation;
- **Moisten skin** before placing the electrodes on the skin;
- Reduce electrode resistance by using the largest electrodes possible for a given anatomical area;
- Optimize electrode placement using the placements recommended in this manual, and;
- If possible, have the person to be treated assist by also voluntarily contracting the stimulated muscles.

Dry skin is highly resistant to the conduction of electric current. Wet (or perspiring) skin has significantly lower resistance. It is true that individuals tend to vary in their normal levels of skin

moisture. It is also best to precede treatments with the application of some form of heat to the areas to be treated. There are a number of means for doing this, from thermal wraps to steam baths or hot whirlpool treatments.

Oily skin also prevents optimal conductivity and should be treated with soap and water prior to start of a treatment session. As a minimum, the skin should be moistened with a sponge or cloth prior to treatment.

D. Application and Handling of Electrodes

When voluntary muscle movements occur, the complex contraction patterns involving several muscles, bones and joints are represented in the cerebral cortex of the brain rather than the movement of each individual muscle. Thus, the brain typically calls a combination of muscles into action; this should be the objective when using neuromuscular stimulation (NMS) as well.

There are several different methods for placing electrodes on the body. The following guidelines are based on using NMS to exercise groups of muscles involuntarily in the same fashion that these muscles contract voluntarily.

Electrode in pairs layout is the recommended electrode placement. Each pair of electrodes is placed on the same side of the body (e.g., one electrode on the right leg vastus lateralis and the other on the same leg on the vastus medialis), both being a part of the quadriceps muscle group.

Split electrode layout is when one pair of electrodes is split on each side of the body (e.g., one electrode on the right quadriceps and one electrode on the left quadriceps). This layout may result in unevenness of contractions which may be difficult to correct.

Bi-polar electrode layout is ideal for large muscles. An electrode is placed at each end of the muscle such that good, well-controlled contractions are obtained.


For all layouts, the electrode placement should correspond to the location of the motor points of the treated muscle(s). For placement guidance, see **Section X**.

THE BOTTOM LINE?

The **Swiss_{STIM} Physio** is a truly valuable adjunct for the medical professional in assuring effective treatments for:

- Prevention or retardation of muscle disuse atrophy
- Relaxation of muscle spasms
- Increasing local blood circulation
- Muscle re-education and strengthening
- Immediate post-surgical stimulation of calf muscles to prevent venous thrombosis
- Maintaining or improving range of motion of extremities

V. Swiss_{STIM} Unit Components



Note: Battery compartment and SMART CARD access located in back of Swiss_{STIM} Physio

VI. Operating Instructions

A. Key Symbols and Functions

S Select Key

Use this KEY to select a program, select a program "phase" and to **START** a selected program. Once a program has started, use this KEY to make further selections in the program.

↑ Multiple Function/Scroll Keys


- Use either key to turn the **Swiss_{STIM} Physio ON**
- When the **Swiss_{STIM} Physio** is ON, use these keys to scroll through the programs
- Once a stimulation program has started, use these keys to **increase** the stimulation *Intensity* (one key for each channel).

↓ Multiple Function/Scroll Keys


- When the **Swiss_{STIM} Physio** is ON, use these keys to scroll through the programs
- Once a stimulation program has started, use these keys to **decrease** the stimulation *Intensity* (one key for each channel)
- After stimulation program is completed use either key to turn the **Swiss_{STIM} Physio OFF**

B. Display (LCD Screen) Indications

Swiss_{STIM} Physio START SCREEN



Example Program Screen



△▽ MENU NAVIGATION KEY SCROLL INDICATORS

P6 PROGRAM NUMBER INDICATOR

Program 1 PROGRAM NAME INDICATOR

20 min stim PROGRAM DURATION INDICATOR (minutes)

Phase 03 PROGRAM PHASE NUMBER INDICATOR



20:00 **SELECTED PROGRAM PHASE DURATION INDICATORS (minutes)**



11 min

WARM-UP PHASE AND DURATION INDICATOR

Indicate that warm-up phase is selected (dark background) and lasts eleven minutes



05 min

SPECIFIC WORK PHASE INDICATOR AND DURATION

Indicate the specific work phase and five minutes duration



04 min

FINAL COOL-DOWN PHASE INDICATOR AND DURATION

Indicate the winding-down phase and four minutes duration



BATTERY CHARGE INDICATOR

Example



80%




93%



STIMULATION INTENSITY INDICATORS ON EACH CHANNEL

Indicate stimulation *Intensity* level (independently set on the two channels)

Example



END OF PROGRAM INDICATOR

The session is over, 4 beeps can be heard. Switch the unit OFF by pressing one of the two keys with a downward pointing arrow (↓)



SMART CARD ABSENT INDICATOR

The smart card on which the programs are recorded is not in place, or is not properly inserted. Switch unit OFF, open the back battery slot and insert the smart card or check that card is correctly positioned.



SMART CARD INVALID INDICATOR

The smart card inserted in the back slot is not valid. Switch unit OFF, open the back battery slot and insert a valid smart card.

INFORMATION ON FUNCTIONS DURING TREATMENT



PAUSE PAUSE: Pause mode during treatment



STOP: Enables user to stop the program before pre-set time



RESUME: Enables user to quit PAUSE mode and resume treatment from where it was stopped.



CHANGE PHASE OPTION: Enables the user to return to the phase menu and repeat or lengthen one particular phase or skip to the following phase.

C. Inserting New Batteries

Slide back cover until battery compartment is open. Insert new battery; use ONLY a 9 volt alkaline or lithium (high quality) battery and ensure proper polarity (+/-).

D. Inserting or Replacing "Smart Card"

The "Smart Card" contains the specific program information for your **Swiss_{STIM} Physio** (your stimulator comes with this card installed). If the card needs to be installed, do so ONLY when your **Swiss_{STIM} Physio** is OFF. Open battery compartment and insert the "Smart Card" into its slot, making sure that the card is properly positioned (LABEL on the card facing up).

E. Turning the Swiss_{STIM} Physio ON

Press one of the scroll keys marked  to turn the **Swiss_{STIM} Physio** ON. This activates the LCD screen and sound but not the *Intensity* function. This function is only activated

after a program selection has been made. If the **Swiss_{STIM} Physio** logo does not appear or if the display indicates low battery strength, the **Swiss_{STIM} Physio** will not function properly. If the message "SWITCH OFF NOW AND INSERT CARD" appears, turn the **Swiss_{STIM} Physio** OFF by pressing one of the ↓ scroll keys. Insert "Smart Card". If the message "INVALID CARD CHANGE CARD" appears, an invalid card has been inserted.

F. Selecting a Stimulation Program and Duration

Make sure the **Swiss_{STIM} Physio** is ON. The LCD display will indicate the number of programs. Scan the programs by using the scroll keys and select desired program and duration. Once a program is selected by using the **S Select** KEY, the LCD display will show specific information about the program (name, total program duration and program phases/phase length). The program will start with the first program phase (unless a specific phase is selected) and will complete ALL phases of the program unless the user selects only a specific program phase.

G. Selecting a "Phase" of a Stimulation Program

Selecting only a program phase is possible. Once a program has been selected, icons for the program phases will appear on the LCD display. Use the ↑ scroll keys to highlight the desired phase and press the **S Select** KEY to select the phase.

H. Automatic Program OFF Function

If no scroll key is pressed within 60 seconds, the **Swiss_{STIM} Physio** will automatically switch OFF to conserve battery life.

I. Selecting Program Intensity

Once a program has been selected, the *Intensity* controls on the **Swiss_{STIM} Physio** are activated. There are 2 scroll keys (a ↑ and a ↓) for each (right and left) channel of the **Swiss_{STIM} Physio**. Use the ↑ scroll key for each channel to **increase** Intensity and the ↓ scroll key to **decrease** the Intensity level. The *Intensity* scale is a 100 step scale, with 100 % corresponding to the maximum *Intensity* level.

J. Using "Pause" Function

Any **Swiss_{STIM} Physio** program can be paused by pressing the **S Select** KEY once. When this is done, the following choices can be made:



STOP – Stops program



RESUME – Program resumes from point it was stopped



CHANGE PHASE – Returns user to program phase menu and allows user to repeat a phase, skip to next phase or repeat the entire program

K. Switching the Swiss_{STIM} Physio OFF During Treatment.

Press both ↓ scroll keys until the *Intensity* level falls to 0%; press one of the scroll keys again and the unit will beep and switch OFF.

L. Connections

Please refer, as needed, to the control schematic of the **Swiss_{STIM} Physio** (page 9)

1. **Ensure that the unit is OFF.**
2. Connect the skin pads (electrodes) to each output cable.
3. Position the skin pads (electrodes) on the motor points of the muscles to be treated (refer to page 23 in this manual for proper positions).
4. Connect the output cables to the **output socket** of the **Swiss_{STIM} Physio** unit. The two (2) cables can be connected to either output socket on the **Swiss_{STIM} Physio**; the **Swiss_{STIM} Physio** unit will function properly

M. Automatic Switch-OFF

The **Swiss_{STIM} Physio** will automatically switch OFF three (3) minutes after completion of a program. To turn the unit OFF immediately, press one of the ↓ scroll keys.

N. Troubleshooting

If the **Swiss_{STIM} Physio** unit is not working, please check the following:

1. Is the battery correctly inserted?
2. Are the cable connectors properly inserted into the **Swiss_{STIM} Physio** unit?
3. Are the skin pads (electrodes) connected to the cables?
4. Are the skin pads (electrodes) adhering to the skin? If not, wet the pad surfaces sparingly with water.

With incorrect, non-symmetrical, stimulation feeling please check the following:

1. Electrodes may be dry or greasy
2. Electrodes not positioned correctly
3. Weak battery
4. Non-symmetrical positioning of the electrodes (See Section X.)
5. Intensity on the two channels are set differently

VII. PROGRAM DESCRIPTIONS

The **Swiss_{STIM} Physio** has four basic muscle stimulation programs; each program is offered in three treatment times for a total of 12 discrete program options. The four basic programs address specific stimulation application for skeletal muscles (see Guide on page 26).

Always stimulate muscles isometrically; make sure that the limb whose muscle(s) is being stimulated is firmly secured to prevent the movement of the limb (resulting from muscle contraction).

Medical professionals should select the appropriate treatment time: 20, 30 or 40 minutes. If longer treatment times are desired, for example if home treatments are prescribed for a patient, the prescribed programs may be repeated by the user or multiple programs used.

Using complete program stimulation cycle (3 phases) is recommended.

Swiss_{STIM} Physio treatment modalities are:

A. RELAXATION OF MUSCLE SPASM

If a muscle spasm is present in any part of the body, the spasm can be relieved by immediate application of stimulation to the affected muscle(s). A spasm, which is usually painful, is inhibited after stimulation and the relief will last for some time. It is normal to experience a return of pain. It may even be reported as an aggravation of pain by some patients. This is only relative and temporary, and the spasm free periods following treatment will become longer and longer until the pain disappears.

To relax the spasms, place at least two electrodes on each muscle in such way that the direction of stimulation (a straight imaginary line between the electrodes) is aligned with the longitudinal axis of the muscle.

Program(s)

Program 1 is best suited; however, Program 4 can be used, if preferred by the patient.

Treatment frequency

No more than every 2 hours, for as long as needed

Stimulation intensity

During the first 2 treatments at light intensity, well below maximum. Thereafter, increase the intensity for subsequent treatments in accordance with patient tolerance.

B. PREVENTION OR RETARDATION OF MUSCLE ATROPHY

When immobilized for medical reasons, all muscles are subject to atrophy (wasting). To remedy this, stimulation of muscles is recommended when normal physical exercise is impossible.

Program(s)

Program 2 only for at least first 5 treatments; follow with Program 2 in combination with Program 3 for more intense workout.

Treatment frequency

Every 24 hours but not less than every 48 hours, for as long as immobilization lasts.

Stimulation intensity

Initially low, progressively increasing to reach the maximum patient tolerance level.

C. PARTIAL ATROPHY AND/OR SPASM OF PARASPINAL MUSCLES

The **Swiss_{STIM} Physio** is especially effective for relieving spasm and accompanying discomfort in the spinal region. The stimulating electrodes should be placed symmetrically on both sides of the spine (see electrode placement photos) at the level of maximum discomfort or in positions indicated by the prescribing physician. The following treatment procedure is recommended:

Program(s)

Program 1 only for the first 5 treatments at maximum tolerable intensity in order to relieve muscle spasm, followed by Program 2 for subsequent treatments.

Treatment frequency

every 24 hours, for 7 to 10 days

Stimulation intensity

initially low, progressively increasing to reach the maximum patient tolerance level

Normally, spasms, which are usually painful, are inhibited right away following the treatment and the relief lasts for some time. It is normal to experience return of pain after some time following each treatment. It may even be reported as an aggravation of pain by some patients. This is only relative and temporary, and the spasm free periods following treatment will become longer and longer until the discomfort disappears. Stimulation in such cases has not only spasm relieving action but also it may reeducate paravertebral muscles to provide better support for the spine, thus preventing the recurrence of back problems.

D. INCREASING LOCAL BLOOD CIRCULATION

Local blood flow increases with stimulation and reaches a peak within the first 15-20 minutes from the onset of stimulation. Therefore, short stimulation treatments are recommended.

Program

Program 1 is recommended, although blood flow increase happens as a beneficial side effect when using any stimulation program.

Treatment frequency

every 24 hours, or as needed

Stimulation intensity

initially low, progressively increasing to reach the maximum patient tolerance level.

E. MUSCLE RE-EDUCATION

Muscle stimulation may be applied to patients before, during or after physiotherapy re-education. Stimulation may also be applied during and/or after immobilization following orthopedic surgery involving the long bones and joints. After deciding which parts of the body are to be treated, and, upon selection of the appropriate number of electrodes and their location, use the following simplified application method:

Program(s)

If muscle spasm is present, start with Program 1 followed by Program 2. Otherwise, use only Program 2.

Treatment frequency

Daily or on alternate days

Number of treatments : Normally 20 treatments, but can be increased to 40 when necessary to reach the treatment goals or recovery of normal or desired functions.

Stimulation intensity

Initially low, progressively increased after two or three treatments, reaching maximum setting, when possible, by the eighth treatment session.

For most effective reeducation of muscles and joints, the following factors are important:

- Try at all times to use the optimal stimulation intensity, i.e., the strongest contractions should be achieved within the limits of patient comfort and tolerance with no pain.
- Stimulation intensity has no absolute value and may vary from day to day and also during a treatment. Thus, when the skin and muscles are warmed (using thermal wraps, for example), or a series of massages (or a steam bath) have been performed before a stimulation treatment, a slightly lower intensity will suffice to give the desired effect.
- It is possible to treat several areas at once on one patient, allowing considerable saving of time.
- It is possible to mount the stimulator electrodes on the skin under a cast and thus prevent muscle atrophy even on a limb which is fully immobilized, as for instance in patients in a pelvopaedic cast.

F. CALF MUSCLE STIMULATION FOR THE PREVENTION OF VENOUS THROMBOSIS BY IMMEDIATE POSTSURGICAL STIMULATION

Muscle stimulation is effective in reducing the evidence of, and preventing, the symptoms of venous thrombosis. After electrode placement on motor points of calf muscles, the following treatment method is recommended:

Program

Program 2

Treatment frequency

As prescribed by physician or up to 3 x per day, for as long as the patient is immobilized

Stimulation intensity

Initially low, progressively increasing up to the maximum patient tolerance level.

G. MAINTAINING OR INCREASING RANGE OF MOTION OF EXTREMITIES

Electrical muscle stimulation may be applied to patients as an adjunct to physiotherapy procedures.

After deciding which extremities are to be treated and upon selection of the appropriate number of electrodes and their location, use the following simplified application method:

Program(s)

If muscle spasm is present, start with Program 1 and follow with Program 3. Otherwise, use only Program 3.

Treatment frequency

Once per day, daily or on alternate days

Number of treatments: Normally 20 treatments but can be increased when necessary to reach the treatment goals or recovery of normal or desired function.

Stimulation intensity

Initially low, progressively increased after two or three treatments, reaching maximum setting, when possible, at the eighth treatment session.

For most effective maintenance of range of motion of extremities, the following factors are important:

- At all times try to use the optimal stimulation intensity; i.e., the strongest contractions within the limits of patient comfort and tolerance, with no pain.
- Stimulation intensity has no absolute value and may vary from day to day and also during a treatment. Thus, when the skin and muscles are warmed (using thermal wraps, for example), or a series of massages (or a steam bath) have been performed before the stimulation treatment, a slightly lower intensity will suffice to give the desired effect.
- It is possible to treat several areas at once on one patient, which permits considerable time savings.

The timing sequences and frequencies for the **Swiss_{STIM} Physio** programs are shown in this table:

		Program 1	Program 2	Program 3	Program 4
Ramp just prior to contraction	sec		3	3	3
Tetanic Contraction-Stimulus ON	sec		6	12	9
Ramp plus contraction	sec		9	15	12
Ramp after contraction	sec		-	-	-
Relaxation after contraction Stimulus OFF	sec		21	39	48
Stimulus frequency, maximum with tetanic contractions	Hz Max	120	87	62	89
Max Impulse duration	µs	300	220	250	180
Total Cycle	sec		30	54	60
Program time (selectable)	min	20/30/40	20/30/40	20/30/40	20/30/40

VIII. Useful Information

A. *Inserting or Replacing Batteries*

The battery charge indicator will indicate when the battery weakens and requires replacement. If the battery is not replaced, the **Swiss_{STIM} Physio** will cease operating within one (1) hour.

Remove (slide) the cover of the unit. Remove old battery and insert a fresh 9-volt lithium, alkaline or rechargeable NiMH battery. Ensure proper battery polarity; the **Swiss_{STIM} Physio** will not operate if polarity is reversed.

B. *Warnings Concerning Battery Handling*

Always read and follow the specific instructions provided by battery manufacturers. Note the following:

- Ensure that battery polarity is correct.
- Do not expose batteries to temperatures exceeding manufacturer's specifications.
- Do not store and/or ship this unit with batteries inserted.
- Do not attempt to recharge alkaline or lithium batteries.
- Do not dispose of any battery in fire.
- Note that batteries may present burn or fire hazard if short-circuited.
- Improper battery handling may result in explosion, leakage or flames.

C. *Handling/Cleaning the Swiss_{STIM} Physio Unit*

Use soft brush or soft cloth to clean unit case; do not use liquid cleansers. Use same procedure with electrical leads. The **Swiss_{STIM} Physio** is designed to be maintenance free.

D. *Storage Conditions*

The **Swiss_{STIM} Physio** may be stored for prolonged periods with no degradation. REMOVE the battery when the unit is not used on a regular basis.

E. *Battery Disposal*

Always dispose of batteries in accordance with battery manufacturer instructions.

F. *Warranty*

Free replacement, for two (2) years, except for accessories and shipping charges. Free replacement will be made if defect is in manufacture of unit; free replacement does not apply when damages are related to improper use or abuse of unit.

IX. Swiss_{STIM} Physio Technical and Safety Data

A. Technical Characteristics

Stimulation Channels

Two, independent, isolated channels. Separate *Intensity* control levels for each channel.

Adjustments

5 keys control all stimulation functions and *Intensity* levels.

Display

The graphical LCD 100 X 30 pixels (46x16 mm) signals all stimulation functions, stimulation intensity per channel, length of program, battery charge status etc.

Output

Current waveform (during training): Asymmetrical biphasic with fast rise and zero net current.
Voltage waveform (open circuit): Low voltage, rectangular, compensated monophasic impulse.
Peak open circuit voltage during each impulse: 45 Vp \pm 10%.

Maximum Output at 500 Ω Load

RMS voltage: \leq 5 V (volts) \pm 10%.

Peak open circuit voltage during treatment: 90mA (milliamperes) \pm 10%.

Power Supply

One 9 Volt lithium, alkaline or NiMH rechargeable battery. Battery charge lasts approximately 20 hours, according to stimulation intensity. The output and program timing parameters are stable (\pm 2%) throughout the life of the battery.

Standard Accessories

Four skin pads (nonpolar, self-adhesive, reusable electrodes)

One alkaline 9 Volt battery

2 cables

Carrying case and User Manual.

Size

Unit size is approximately 12.5*8.1*2.6 cm. Weight is 200 grams or approximately 7 oz.

Cleaning/Maintenance

Use soft brush or soft cloth to clean unit case; do not use liquid cleansers. Use same procedure with electrical leads. The **Swiss_{STIM} Physio** is designed to be maintenance free.

B. Safety

The **Swiss_{STIM} Physio** meets the requirements of ANSI/AAMI NS4 1985 American Standard and the IEC-601-2-10 European Safety Standard.

Specific Safety Features and Measures to Prevent Misuse

The **Swiss_{STIM} Physio** is inherently designed so that improper or accidental application of this unit **cannot** produce cardiac rhythm disturbances. ***This does not include individuals with implanted cardiac pacemakers!***

This inherent safety is due to the minimal electric charge of all stimulating impulses in the **Swiss_{STIM} Physio**. In all programs of stimulation, the typical impulse is below 10 micro coulombs and the maximum possible impulse is 20 micro coulombs, zero net current. Moreover, this maximum charge of 20 micro coulombs is significantly below the cardiac rhythm disturbance safety margin of 25 micro coulombs per pulse, the standard as established by the Association for Advancement of Medical Instrumentation (AAMI) in AAMI/ANSI standard NS4-1985 and acknowledged by the U.S. Food and Drug Administration (FDA).


Specific safety features in the **Swiss_{STIM} Physio** include:

1. Impossible for users to modify the embedded programs; users can only modify the intensity of stimulation
2. All programs begin with minimal electrical *Intensity*; the user must increase the *Intensity* to the desired training level.
3. Maximum possible electrical impulse is 20 micro coulombs.
4. Cable design prevent the possibility of accidental connection to a power source, such as an AC power outlet.
5. Automatic control of stimulation current density precludes excessive current density at the electrode-skin interface and ensures skin safety.


X. Electrode Placement

Use oval 3" by 5" electrodes for large areas (e.g., leg muscles) and round 3" electrodes for smaller areas such as forearm muscles. Do not use smaller electrodes than recommended above. **Leads for Channel 1 are depicted in RED and leads for Channel 2 in BLUE.** Where only one set of electrodes is shown, Channel 1 is depicted but Channel 2 may be used instead. Channels may be reversed, if desired, from that indicated in the following photos.

QUADRICEPS




On both legs




On one leg

TIBIA




Channel 2 electrodes may be attached to other leg in order to treat both legs

CALVES




THIGHS




Channel 2 electrodes may be attached to other thighs in order to treat both thighs


FEMORAL BICEPS




FEMORAL BICEPS AND CALVES




FOREARM




BICEPS




DELTOID




TRICEPS




ABDOMINALS




UPPER TRAPEZIUS




OR




LOWER TRAPEZIUS




LUMBAR VERTEBRAE



OR



NECK BACK



XI. Program Reference Guide

Indication to be Treated	Program 1 (P1)	Program 2 (P2)	Program 3 (P3)	Program (P4)
RELAXATION OF MUSCLE SPASM	X			X If preferred by patient
PREVENTION OR RETARDATION OF MUSCLE ATROPHY		X Use P2 for first 5 treatments	X Use in <i>combination</i> with P2 after 5 treatments	
PARTIAL ATROPHY AND/OR SPASM OF PARASPINAL MUSCLES	X Use P1 for at least first 5 treatments	X Use P2 after first 5 treatments		
INCREASING LOCAL BLOOD CIRCULATION	X			
MUSCLE RE-EDUCATION	X Use only if spasm is present, then continue with P2	X		
CALF MUSCLE STIMULATION FOR THE PREVENTION OF VENOUS THROMBOSIS BY IMMEDIATE POSTSURGICAL STIMULATION		X		
MAINTAINING OR INCREASING RANGE OF MOTION OF EXTREMITIES	X Only with spasm present, continue with P3		X	